# Digital Refractometer

for Sodium Chloride Measurement Throughout the Food Industry

- Ideal for the analysis of:
  - Salad dressings, cheeses, condiments, pickles, canned foods, jarred foods, milk, juices, energy drinks, soups, brines and whey
- High accuracy measurements in g/100 g, g/100 mL, specific gravity and °Baume
- Dual-level LCD
  - The dual-level LCD displays measurement and temperature readings simultaneously
- ATC
  - Automatic Temperature
    Compensation
- Easy measurement
  - Place a few drops of the sample in the well and press the READ key
- BEPS
  - Alerts the user of low battery power that could adversely affect readings.
- IP65 water protection
  - Built to perform under harsh laboratory and field conditions
- Quick, accurate results
  Readings are displayed in approximately 1.5 seconds.
- Single-point calibration
  - Calibrate with distilled or deionized water
- Small sample size
  - Sample size can be as small as 2 metric drops
- Automatic shut-off
  - After three minutes of non-use
- Stainless steel sample well
- Easy to clean and corrosion resistant
- ABS thermoplastic casing



### Ideal for the Food Industry

Hanna offers the HI96821 digital sodium chloride refractometer to meet the requirements of the food industry. This optical instrument employs the measurement of the refractive index to determine sodium chloride concentration in aqueous solutions used in food preparation. It is not intended for seawater salinity measurements.

#### **Refractive Index**

The measurement of refractive index is simple and quick and provides the user an accepted method for sodium chloride analysis. Samples are measured after a simple user calibration with deionized or distilled water. Within seconds the instrument measures the refractive index of the solution, apply the necessary calculations and display the results in the selected unit. The digital refractometer eliminates the uncertainty associated with mechanical refractometers and is portable for measurements where you need them.

#### Features

The instrument utilizes internationally recognized references for unit conversion and temperature compensation. It can display the measurement of NaCl concentration 4 different ways: g/100 g, g/100 mL, Specific Gravity, and °Baumé.

Temperature (in °C or °F) is displayed simultaneously with the measurement on the large dual level display along with icons for Low Power and other helpful message codes.

13

nts | www.hannainst.com

## Easy to Operate

#### Startup Screens

When the HI96821 is turned on, all of the LCD segments will be displayed followed by the percentage of battery life remaining.

#### Calibration

Perform a quick and easy calibration after startup:

- 1. Using a pipette, completely cover the prism in the sample well with distilled or deionized water.
- 2. Press the ZERO key.

#### Unit Selection

Just press the RANGE key to cycle through the HI96821's units of measurement (g/100 g, g/100 mL, Specific Gravity and °Baumé).

#### Measurement

Achieve fast, accurate results:

- 1. Using a plastic pipette, place sample onto the prism surface until the well is full.
- 2. Press the READ key and the results are displayed in the selected units.

#### Making a Standard Sodium Chloride Solution

To make a standard NaCl solution (g/100 g), follow the procedure below:

- Place a container (such as a glass vial or dropper bottle that has a cover) on an analytical balance.
- Tare the balance.

Information

- To make an X NaCl solution weigh out X grams of high purity dried Sodium Chloride (CAS #: 7647-14-5: MW 58.44) directly into the container.
- Add distilled or deionized water to the container so the total weight of the solution is 100 g.

Example with g/100 g NaCl:	
g/100 g NaCl	10
g NaCl	10.000
g Water	90.000
g Total	100.000

#### HI96821 Specifications Range 0 to 28 Resolution 0.1 g/100 g Accuracy ±0.2 (@25°C/77°F) 0 to 34 Range Resolution 0.1 g/100 mL Accuracy ±0.2 (@25°C/77°F) 1.000 to 1.216 Range 0.001 Resolution Specific Gravity (S.G.) Accuracy +0.002 (@25°C/77°F) Range 0 to 26 Resolution 0.1 °Baumé Accuracy ±0.2 (@25°C/77°F) 0 to 80°C (32 to 176°F) Range Resolution 0.1°C (0.1°F) Temperature Accuracy ±0.3°C(±0.5°F) (@25°C/77°F) Temperature automatic between 10 and 40°C (50 to 104°F) Compensation Measurement Time approximately 1.5 seconds Minimum Sample 100 µL (to cover prism totally) Volume Additional Light Source yellow LED Specifications Sample Cell stainless steel ring and flint glass prism Auto-off after three minutes of non-use Enclosure Rating IP65 Battery Type / Battery Life 9V / approximately 5000 readings Dimensions / Weight 192 x 102 x 67 mm (7.6 x 4.01 x 2.6") / 420 g (14.8 oz.) Ordering

HI96821 is supplied with battery and instruction manual.

12

